



# BIOM9660 Assignment

## Software Design and more....

Never Stand Still

Faculty of Engineering

Graduate School of Biomedical Engineering

Greg Watkins – August 2016

## ▪ **Software Design & Submission**

- Frequency/Time Matrix
- FFT as a band pass filter
- Amplitude Mapping



Dr B.Swanson – Cochlear Implant Sound Processing model

# Software – Files Provided

- Three MatLab files will be provided
  1. **cochlearProc.m** - runs your Cochlear Implant simulation
  2. **classCochlear.m** - skeleton of functions you need to write
  3. **classCochlearSupport.m** - “helper” functions/data
- These files should be copied to a single folder on your PC
- The only file you may modify is #2 – more later

# classCochlear.m – *Insert your code (only) here!*

- This file provides the declarations for 6 functions that you must write

- The declarations (names and parameters) must not be changed
- Comments provide a brief description of what the function should do. The details are up to you.

- For example

- `function result = getWav(obj, name)`  
*% read the wav file <name> and resample if required so*  
*% that the sampling frequency is fSample +/- fTolerance*

getWav

getFTM

process

applyDR

plotSignal

plotElectrodiagram

# Software –Testing

- To test your simulation in MatLab type....

*cochlearProc( soundFileName, processingType)*

Where

soundFileName is the name of a sound file eg 'hello.wav',  
processingType is 1 or 2 or 3 (see assignment spec)

- This function will call the functions you have written and save the results to disk. Eg output files for “Speak” are:-
  - Speak.jpg (electrodogram)
  - Speak\_speech.jpg (time domain speech)
  - Speak.csv (stimulation file)
  - Speak.wav (vocoded sound file)

# Other data and functions

- `classCochlearSupport.m` provides helper functions and data
- For example
  - `fSample = 16000;` % required sampling frequency
  - `fTolerance = 0.1;` % allowed variation from `fSample`
  - `tSample = 0.002;` % sample time for cochlear implant processing
  - `frameOverlap = 3/4;` % required overlap when using Hann window
  - `numElectrodes = 16;` % number of electrodes in cochlear implant
  - `numFormants = 2;` % formants used by the F0F1F2 strategy
  - `dynamicRange = 10;` % dB
  - `maxOutput = 1024;`
  - `function result = procName(obj, procType)`  
% return string containing the name of the process type <procType>
- Access data using `obj.<name>` eg `obj.numElectrodes`

# Software Submission

- Only your file **classCochlear.m** is submitted via Moodle
- When you submit classCochlear.m you must rename it to 'class<znumber>' (eg classz1234567)
  - Moodle will ask for two versions 'class<znumber>.m' and 'class<znumber>.txt' (eg classz1234567.m and classz1234567.txt)
  - If your files have other names they will not be marked.
- cochlearProc.m and classCochlearSupport.m are not submitted.
  - Any changes will not be used when your code is marked.

- Software Design & Submission
- **Frequency/Time Matrix**
- FFT as a band pass filter
- Amplitude Mapping




Dr B.Swanson – Cochlear Implant Sound Processing model



# Frequency/Time matrix (FTM)

- Every 2mS your simulation will take a sample (frame or epoch) of the input signal and pass it through  $n$  parallel band pass filters (BPF)
- This provides a matrix as shown at right called the “Frequency/time matrix” or FTM
- The FTM is the key data structure in your simulation




- Software Design & Submission
- Frequency/Time Matrix
- **FFT as a band pass filter**
- Amplitude Mapping



Dr B.Swanson – Cochlear Implant Sound Processing model

# FFT – A possible BPF implementation


- As an example... take a sample of a 1kHz Square wave
- Nyquist theorem:  $F_s \geq 2BW$ . For  $F_s=16\text{kHz}$  max bandwidth is 8kHz.




Sampled waveform  
split into “frames”  
and FFT applied.

MatLab FFT

Use abs() function to  
obtain magnitude of  
bin data




# BPF/FFT



In cochlear implants FFT used to implement  $n \times$  bandpass filters

Need to choose number of bins appropriately


The FFT is a simple way to implement the required Band-pass filters



- The example shown has a wave form that fits perfectly in the sample frame. Very neat plots result.
- In real-life and in particular with speech this is not always true....

# Hann Window


- Real-life signals, especially speech, don't fit neatly in sample window.
- Windowing and overlapped sampling are applied



Plot at left shows a wave form that does not fit neatly in the sample window.

Hann window (red) is applied to remove transients

Plot at right shows  $n \times$  overlapping, samples with hann window applied as required in the assignment



- Software Design & Submission
- Frequency/Time Matrix
- FFT as a band pass filter
- **Amplitude Mapping**



Dr B.Swanson – Cochlear Implant Sound Processing model

# Amplitude Mapping



- Referring to Dr Swanson's processing model the “Amplitude Mapping” block contains a number of functions including:-
  1. Ensuring that levels are limited within specific bounds (eg 0 to 1)
  2. Mapping the levels to an individuals T&C levels
- The assignment “stimulus amplitude range” requirement is similar to 1) + 2) above.

# Amplitude Mapping – T & C levels

- The assignment specifies
  - maximum or “C” level of 1024 and
  - dynamic range (difference between C and T levels) of 10dB.
- Rearranging the equation from the assignment instructions gives

$$A_{\min} = \frac{A_{\max}}{\sqrt{10}} = 324 = "T"$$


where  $A_{\max}$  is the “C” level and  $A_{\min}$  is the “T” level

- This tells us that the amplitude information presented to the Vocoder or plotted in the electrodogram must always be in the range 324 to 1024.




# Amplitude Mapping – Scaling

- Let's assume that the amplitude of the information output from the processing strategy in the “Sampling & Selection” block (eg ACE) is in the range 0 to 1.
- In the simplest approach
  - An input amplitude of 0 is mapped to the “T” level
  - An input amplitude of 1 is mapped to the “C” level
  - Amplitudes in between are scaled between the “T” and “C” levels



# Amplitude Mapping - Scaling

- As part of your assignment you might investigate whether more sophisticated mappings are more intelligible.
- The mapping at right maps any low amplitude input (eg  $< 0.2$ ) to 0 and scales inputs  $> 0.2$  between the T and C levels
- This approach is used in Cochlear Implants. You might like to consider the advantages of this approach





**Never Stand Still**